

# Surgical Face Mask



- CE
- CFDA
- COVID-19 Medical Device

Authorization for  
Importation or Sale in  
Canada

- White List

**Material** Three-layer Mask  
 Inner layer: Hygroscopic nonwoven fabric  
 Protective layer: 25gms Melt-blown fabric  
 Outer later: Waterproof nonwoven fabric

**Size** 17.5\*9.5cm

**Standard** EN 14683:2019

Model Name	Type I (Ear loop)	Type II (Ear loop)	Type IIR (Ear loop)	Type IIR (Tie on)
Cat. No.	O06B010110	O06B020110	O06B020210	O06B020310
Bacterial filtration efficiency	>95%	>98%	>98%	>98%
Differential pressure	<40Pa/cm <sup>2</sup>	<40Pa/cm <sup>2</sup>	<60Pa/cm <sup>2</sup>	<60Pa/cm <sup>2</sup>
Splash resistance pressure	/	/	≥16kPa (120mmHg)	≥16kPa (120mmHg)

## WELL LEAD MEDICAL CO., LTD.

Address: 47 Guomao Avenue South, Panyu, Guangzhou, China

Tel: +86-20-84758878 Fax: +86-20-84758224 Post Code: 511434

Marketing & Trading Office

Add: RM 805, No 55, 21 Long, Pusan Road, Shanghai, China

Tel: +86-21-58397623 Fax: +86-21-58397615 Post code: 200125

E-mail: info@welllead.com.cn <http://www.welllead.com.cn>

**wellead**®

# 医 | 用 | 口 | 罩



材 料：三层面罩                      内 层：无纺布  
尺 寸：17.5\*9.5cm                    中间层：熔喷布  
执行标准：EN 14683:2019          外 层：防水无纺布

型号名称	Type I (挂耳式)	Type II (挂耳式)	Type IIR (挂耳式)	Type IIR (绑带式)
目录号	O06B010110	O06B020110	O06B020210	O06B020310
细菌过滤效率BFE	>95%	>98%	>98%	>98%
两侧压力差	<40Pa/cm <sup>2</sup>	<40Pa/cm <sup>2</sup>	<60Pa/cm <sup>2</sup>	<60Pa/cm <sup>2</sup>
抗喷溅压力	/	/	≥16kPa (120mmHg)	≥16kPa (120mmHg)

## 广州维力医疗器械股份有限公司

地址: 广东省广州市番禺区国贸大道南47号  
电话: +86-20-84758878 传真: +86-20-84758224 邮政编码: 511434  
上海维力医疗用品进出口公司  
地址: 上海市浦东新区浦三路21弄55-56号银亿滨江中心壹号楼805室  
电话: +86-21-58397623 传真: +86-21-58397615 邮政编码: 200125  
邮箱: info@welllead.com.cn 网址: http://www.welllead.com.cn

welllead®