

Surgical Face Mask



- CE
- CFDA
- COVID-19 Medical Device

Authorization for
Importation or Sale in
Canada

- White List

Material Three-layer Mask
 Inner layer: Hygroscopic nonwoven fabric
 Protective layer: 25gms Melt-blown fabric
 Outer later: Waterproof nonwoven fabric

Size 17.5*9.5cm

Standard EN 14683:2019

| Model Name | Type I (Ear loop) | Type II (Ear loop) | Type IIR (Ear loop) | Type IIR (Tie on) |
|---------------------------------|-----------------------|-----------------------|-----------------------|-----------------------|
| Cat. No. | O06B010110 | O06B020110 | O06B020210 | O06B020310 |
| Bacterial filtration efficiency | >95% | >98% | >98% | >98% |
| Differential pressure | <40Pa/cm ² | <40Pa/cm ² | <60Pa/cm ² | <60Pa/cm ² |
| Splash resistance pressure | / | / | ≥16kPa (120mmHg) | ≥16kPa (120mmHg) |

WELL LEAD MEDICAL CO., LTD.

Address: 47 Guomao Avenue South, Panyu, Guangzhou, China

Tel: +86-20-84758878 Fax: +86-20-84758224 Post Code: 511434

Marketing & Trading Office

Add: RM 805, No 55, 21 Long, Pusan Road, Shanghai, China

Tel: +86-21-58397623 Fax: +86-21-58397615 Post code: 200125

E-mail: info@welllead.com.cn http://www.welllead.com.cn

wellead[®]

医 | 用 | 口 | 罩



材 料：三层面罩 内 层：无纺布
尺 寸：17.5*9.5cm 中间层：熔喷布
执行标准：EN 14683:2019 外 层：防水无纺布

| 型号名称 | Type I (挂耳式) | Type II (挂耳式) | Type IIR (挂耳式) | Type IIR (绑带式) |
|-----------|-----------------------|-----------------------|-----------------------|-----------------------|
| 目录号 | O06B010110 | O06B020110 | O06B020210 | O06B020310 |
| 细菌过滤效率BFE | >95% | >98% | >98% | >98% |
| 两侧压力差 | <40Pa/cm ² | <40Pa/cm ² | <60Pa/cm ² | <60Pa/cm ² |
| 抗喷溅压力 | / | / | ≥16kPa (120mmHg) | ≥16kPa (120mmHg) |

广州维力医疗器械股份有限公司

地址: 广东省广州市番禺区国贸大道南47号
电话: +86-20-84758878 传真: +86-20-84758224 邮政编码: 511434
上海维力医疗用品进出口公司
地址: 上海市浦东新区浦三路21弄55-56号银亿滨江中心壹号楼805室
电话: +86-21-58397623 传真: +86-21-58397615 邮政编码: 200125
邮箱: info@welllead.com.cn 网址: http://www.welllead.com.cn

wellead®