

● Capnography Sampling Mask

■ Features

- Designed for oxygen supply and sampling of expired CO₂ gas simultaneously
- Offered with head strap and adjustable nose clip
- The star lumen tubing can ensure oxygen flow even if the tube is kinked
- The standard length of tube is 2.1m, and different length is available

Type	Specification	Cat. No. (PVC)	Cat. No. (DEHP-FREE)
Adult standard	2.1m tubing	R01A100110	R01A10011F
Adult elongated	2.1m tubing	R01A100210	R01A10021F
Pediatric standard	2.1m tubing	R01A100310	R01A10031F
Pediatric elongated	2.1m tubing	R01A100410	R01A10041F

