

● Non-Rebreathing Mask

■ Features

- Provide a high concentration of oxygen
- Offered with head strap and adjustable nose clip
- The star lumen tubing can ensure oxygen flow even if the tube is kinked
- The standard length of tube is 2.1m, and different length is available

Type	Specification	Cat. No. (PVC)	Cat. No. (DEHP-FREE)
Adult standard	2.1m tubing	R01A040110	R01A04011F
Adult elongated	2.1m tubing	R01A040210	R01A04021F
Pediatric standard	2.1m tubing	R01A040310	R01A04031F
Pediatric elongated	2.1m tubing	R01A040410	R01A04041F
Infant	2.1m tubing	R01A040510	R01A04051F

