

## Tracheostomy Tube

- Tracheostomy Tube is a hollow tube, with or without a cuff, that is electively inserted directly into the trachea through a surgical incision or with a wire-guided progressive dilatation technique in case of an emergency.
- Lightweight tube made of Medical grade PVC for patient comfort. Full length Radio-opaque line for detection of correct placement. ISO standard connector for universal connection to ventilation equipment Printed neck plate with size information for easy identification.
- Straps provided in the pack for fixation of the tube. Smooth rounded tip of the Obturator reduces trauma during insertion. High volume, low pressure Cuff provides excellent sealing. Rigid blister pack provides maximum protection for the tube.



### ■ Cuffed

Size ID (mm)	Cat. No.	Cat. No. (DEHP-FREE)
3.0	A02A033010	A02A03301F
3.5	A02A033510	A02A03351F
4.0	A02A034010	A02A03401F
4.5	A02A034510	A02A03451F
5.0	A02A035010	A02A03501F
5.5	A02A035510	A02A03551F
6.0	A02A036010	A02A03601F
6.5	A02A016510	A02A01651F
7.0	A02A017010	A02A01701F
7.5	A02A017510	A02A01751F
8.0	A02A018010	A02A01801F
8.5	A02A018510	A02A01851F
9.0	A02A019010	A02A01901F
9.5	A02A019510	A02A01951F
10.0	A02A019610	A02A01961F



### ■ Uncuffed

Size ID (mm)	Cat. No.	Cat. No. (DEHP-FREE)
3.0	A02A043030	A02A04303F
3.5	A02A043530	A02A04353F
4.0	A02A044030	A02A04403F
4.5	A02A044530	A02A04453F
5.0	A02A045030	A02A04503F
5.5	A02A045530	A02A04553F
6.0	A02A046030	A02A04603F
6.5	A02A026530	A02A02653F
7.0	A02A027030	A02A02703F
7.5	A02A027530	A02A02753F
8.0	A02A028030	A02A02803F
8.5	A02A028530	A02A02853F
9.0	A02A029030	A02A02903F
9.5	A02A029530	A02A02953F
10.0	A02A029630	A02A02963F



### ■ Tracheostomy tube with PU Cuff

Size ID (mm)	Cat. No.	Cat. No. (DEHP-FREE)
3.0	A02B023040	A02B02304F
3.5	A02B023540	A02B02354F
4.0	A02B024040	A02B024F4F
4.5	A02B024540	A02B02454F
5.0	A02B025040	A02B02504F
5.5	A02B025540	A02B02554F
6.0	A02B026040	A02B02604F
6.5	A02B016540	A02B01654F
7.0	A02B017040	A02B01704F
7.5	A02B017540	A02B01754F
8.0	A02B018040	A02B01804F
8.5	A02B018540	A02B01854F
9.0	A02B019040	A02B01904F
9.5	A02B019540	A02B01954F
10.0	A02B019640	A02B01964F

### ■ SuctionPlus Tracheostomy Tube

Size ID (mm)	Cat. No.	Cat. No. (DEHP-FREE)
6.0	A02C016010	A02C01601F
6.5	A02C016510	A02C01651F
7.0	A02C017010	A02C01701F
7.5	A02C017510	A02C01751F
8.0	A02C018010	A02C01801F
8.5	A02C018510	A02C01851F
9.0	A02C019010	A02C01901F